


Main specifications

Spatial Resolution		1.2 nm @ 30 kV 3.0 nm @ 1 kV 2.0 nm @ 1 kV with deceleration mode*1 3.0 nm @ 15 kV Variable Pressure mode*2	
Magnification		10~600,000 × (based on 4 " × 5 " picture) 18~1,000,000 × (800×600 pixels on display) 30~1,500,000 × (1,280×960 pixels on display)	
Electron Optics	Electron gun	ZrO/W Schottky emission electron source	
	Acceleration Voltage	0.5~30 kV (0.1 kV step)	
	Landing Voltage	0.1~2.0 kV with deceleration mode*1	
	Maximum Probe Current	> 200 nA	
Detector		Evehart Thonley SE detector (Lower detector)	
Variable Pressure Mode *2		Pressure Range: 10~300 Pa	
Specimen Stage	Control		
	Movement	X	0~100 mm
		Y	0~50 mm
		Z	3~65 mm
		T	-20~90 °
		R	360 °
Specimen Size		up to 200 mm φ maximum 80 mm height	
Monitor *3		23 inch LCD (1,920×1,080)	
Display Mode	Large screen display	1,280×960 pixels	
	Single image display	800×600 pixels	
	Dual image display	640×480 pixels	
	Quad image display	640×480 pixels	
Image Data Saving	Pixel Size	640×480 1,280×960 2,560×1,920 5,120×3,840	
Dimension and Weight	Main Unit	740(W) × 1,000(D) × 1,650(H) mm 550 kg	
	Display	1,100(W) × 1,120(D) × 730(H) mm 290 kg	
	Rotary Pump *3	526(W) × 225(D) × 306(H) mm 28 kg	
	Air Compressor *3	400(W) × 230(D) × 520(H) mm 18 kg	
	Weight	200(W) × 180(D) × 160(H) mm 40 kg	
Optional Detectors	Top detector for high resolution imaging Through-the-Lens*1 Ultra Variable-Pressure Detector (UVD) Retractable five segment Backscatter Electron Detector (PD-BSD) *4 Energy Dispersive X-ray detector (EDS) Wavelength Dispersive X-ray detector (WDS) Electron Backscatter Diffraction Pattern detector (EBSD) CL detector STEM detector		

Automatic functions	Auto calibration function		
	Auto Start		
	Automatic ABCC		
	Automatic AFC		
	Automatic ASF		
	Auto Beam Blanking		
	Operator assist functions	EM Wizard	Two separated mode selectable Assisting observation purpose
		Astigmatism learning	Astigmatism correction training
		Operation guide	Step-by-step operation flow guide
		Application assist	Pop-up help for suggesting cause and solution
Tilt assist		Assisting tilt function for finding EBSD analysis available area	
Stage memory			
Move stage		Acquired image only	
Image navigator			
Condition save files			
NTSC video output			
Optional Softwares	Movie capture		
	Live stereoscopic function		
	Tilt compensation function		
	CD measurement function	Offline CD measurement	
	Zigzag capture		
	Stitch software		
Optional accessories	Cryo stage with sample preparation chamber		
	Heating stage		
	Cooling stage		
	Airlock		
	Plasma cleaner		
	Pico ampere meter		
	Faraday cup		
	Multi specimens holder		

- *1: Top detector is option, combined to deceleration function.
*2: Variable Pressure mode is option.
*3: option
*4: PD-BSD is standard detector in Variable Pressure system.



Notice: For correct operation, follow the instruction manual when using the instrument.

Specifications in this catalog are subject to change with or without notice, as Hitachi High-Technologies Corporation continues to develop the latest technologies and products for our customers.

Copyright (C) Hitachi High-Technologies Corporation 2015 All rights reserved.

Hitachi High-Technologies Corporation

Hitachi High Technologies America, Inc.

Toll free: 800-548-9001 (US&Canada)

www.hitachi-hightech.com

Schottky FE-SEM SU5000

HITACHI
Inspire the Next

SU5000


FIELD EMISSION / VARIABLE PRESSURE
SCANNING ELECTRON MICROSCOPE




High resolution at low accelerating voltage

2.0 nm SE Image resolution at 1 kV

The SU5000 has a redesigned Schottky field emission electron gun, which produces a very fine electron beam, with high brightness and a narrow energy width. Combined with a new, low aberration objective lens, beam deceleration and an in-column detector, the SU5000 can capture high resolution images at voltages as low as 100 V.



Landing Voltage: 1 kV, Secondary Electron (SE) Image
Magnification: 120,000 x, Sample: Zinc oxide powder



Landing Voltage: 500 V, Secondary Electron (SE) Image
Magnification: 70,000 x, Sample: Anodized Aluminum oxide

High probe current with high performance

Advanced optics for >200 nA of probe current and simple switching between high and low vacuum modes

The SU5000's unique design allows for a simple transition between high vacuum and variable pressure modes. There is no need to physically change any pressure limiting apertures inside the chamber, so there is no risk of damaging anything. No aperture change also means no reduction in the imaging field of view or loss of probe current, when going from high vacuum to variable pressure. The SU5000 always provides >200 nA of current, perfect for EDX, WDX, EBSD, and CL.

	Low magnification : 10 x	Probe current: 50 nA Magnification : 50,000 x	Probe current: 200 nA Magnification : 50,000 x
HV High Vacuum < 7 x 10 ⁻⁴ Pa			
LV Low Vacuum 15 Pa			
	No reduction in field of view between HV and LV modes	High probe current images shown at 50 nA and 200 nA for analytical work under both HV and LV modes	

Unique five segment BSE detector

Hybrid concentric+Quad BSE detector distinguishes high angle, low angle and directional BSE signals


Concentric detectors can separate low and high angle BSEs, but have no directional information. 4-quadrant detectors are just the opposite. The SU5000's new hybrid BSED does both, revealing crystal orientation and surface (left) or compositional details (right), depending on the signal orientation to the electron beam axis.



BSED graphic showing the single inner segment and four independent outer segments



Accelerating Voltage: 15 kV, Outer segment/High angle BSE
Magnification: 7,000 x, Sample: Heat resistant steel
This accentuates the crystal orientation of the material




Accelerating Voltage: 15 kV, Inner segment/Low angle BSE
Magnification: 3,000 x, Sample: Heat resistant steel
This accentuates the chemical composition of the material


Charge Suppression Scan Technology

CS Scan for charge reduction and limiting radiation damage for beam sensitive samples

CS Scan changes the way the electron beam moves across the sample, reducing time the primary beam dwells at a location. This helps reduce sample charging and radiation damage on beam sensitive samples.



Accelerating Voltage: 2 kV, SE image
Magnification: 30,000 x
Sample: Aluminum electrolytic capacitor
Slow scan: 32 seconds
Without CS Scan




Accelerating Voltage: 2 kV, SE image
Magnification: 30,000 x
Sample: Aluminum electrolytic capacitor
Slow scan: 32 seconds
With CS Scan


VP Mode Detection Technology

Ultra Variable-pressure Detector (UVD) for high quality, SE-like images in variable pressure mode

The UVD detects photons emitted from collisions between electrons and gas molecules in the sample chamber, creating images with excellent surface detail and topographic information in variable pressure mode. The UVD is efficient at low and high accelerating voltages, across the full pressure range of the microscope, making it a versatile tool for variable pressure imaging.




Accelerating Voltage: 3 kV, Variable Pressure, 30 Pa, UVD
Magnification: 7,000 x, Sample: PTFE



UVD Principle

“EM Wizard”

EM wizard is a completely new, knowledge-based system for SEM imaging that goes beyond basic pre-set conditions and recipes. Self-directed optical alignment, fast and accurate automated image adjustment functions and a simplified user interface redefine ease of use for a field emission SEM.




Any operator can acquire quality images with a few point-and-click selections, telling the SEM the intended purpose for observing the sample. Integrated visual applications assistance, practical guidance, tutorials and training tools inspire the user to learn more.

Novice or expert, the results can now be the same.

Elemental + Surface information
Specimen : Neodymium magnet

Specimen : Ferrite core
Surface information




SEM image	Radar chart	Observation purpose
		Standard observation
		Surface information
		Elemental information
		Elemental + Surface information
		Element analysis + Observation

*above table : BSD fitted configuration. BSD is option.

Intuitive user interface

Standard mode offers simple and quick operation. Novice users are guided, step by step, and learn along the way, thanks to the interactive user guide.



The advanced mode provides full access to all SEM controls the expert requires, with versatile control and display of multiple detector signals. The advanced 3D MultiFinder stage navigation safely carries out complex tilt and rotation stage movements, based on simple click and drag input. 3D MultiFinder also provides precise positioning of the sample's field of view for high tilt EBSD analysis.

Advanced guide


Guides the user step-by-step and educates the operator along the way.



Powerful Automated Alignment


Auto calibration (US patent; #6864493)

Beam axis and astigmatism alignments are automatically calibrated, stored and then adjusted for a variety of SEM conditions, corresponding to different EM applications. It can also restore the SEM back to its "best condition" state whenever necessary.




Operation tutorial

Teaches proper focus and astigmatism correction via simulations.



Camera navigation


Camera navigation lets the operator quickly move the stage to an area of interest, using a low magnification digital image. The image is most often from the integrated navigation camera but can also be a low mag SEM image or imported from any digital camera source.




Large area SEM image*1

Multiple SEM images are automatically collected and saved to disk. Subsequent wide area SEM images are created by "stitching" together the stored images.

Metal cross section
Accelerating Voltage :
15 kV,
Magnification :
700 x





Total 192 images (16 x 12) automatically stored.



Wide area SEM image can be created after stitching 192 images.

Live Stereoscopic Imaging Function

Live stereoscopic imaging enables real-time 3-D SEM viewing, without tilting the specimen. Live stereoscopic images are generated by rapidly alternating the electron beam tilt angle, to yield left and right parallax images. The parallax images are then synchronized and observed directly with colored 3D eyeglasses.



Analytical Versatility

Multi-functional Specimen Chamber


The large SU5000 analytical chamber has 11 ports to accommodate various accessories and supports simultaneous EDX, WDX, and EBSD. The drawer style chamber with door mounted stage provides for safe and easy sample exchanges, with a fast 3 minute pump down time.

EDS^{*1}


Sample: Bonding solder (Cross-section)
Accelerating voltage: 15 kV, Ip: 2 nA
Acquisition time: 5 min




The pictured EDX and EBSD and WDX components are optional




Backscatter electron image



Copper EDS map (Cu-L)



Lead EDS map (Pb-M)




Tin EDS map (Sn-L)

Cryo stage^{*1}


Cryo preparation for SEM is a common form of observation for "beam sensitive" specimens.



The pictured Cryo-SEM system is optional




Cryo stage attachment




Accelerating Voltage: 2 kV, SE image
Magnification: 2,000 x, Sample: Brevipalpus (Mite)
Cryo temperature: -120 °C
Brevipalpus species have different surface contours.
Sample courtesy of USDA, ARS - Beltsville, MD
Drs. Gary Bauchan & Ronald Ochoa


EBSD^{*1}




EBSD requires a sample to be tilted to 70 deg. This high tilt condition causes substantial trapezoidal image distortion of rectangular areas, due to changes in focus, magnification, beam rotation and other factors. Traditional routines, such as dynamic tilt and focus compensation, attempt to correct for these distortions, but close inspection reveals that spatial inaccuracies remain. Hitachi has developed a new solution that corrects for all conditions to maintain the precise shape and dimensions of the original area. For EBSD, this is essential for reliable grain measurements and accurate stitching of multi-field orientation maps.




70 deg tilt
Without tilt compensation
Magnification: 30 x




70 deg tilt image
With tilt compensation
Magnification: 30 x




10 x 8 montage at 70 deg tilt.
Total Area is Approx. 6 mm x 6 mm
Showing precisely corrected trapezoid and accurate stitching.



BC map



IPF-Z map



Multi-field EBSD map of heat resistant steel, 11 x 9 matrix, 2 um step, total area 5 x 3.8 mm.


2mm

2mm


200 um

Cooling stage^{*1}


Peltier cooling of specimens (0 to -50 °C) can be utilized to slow the evaporation process under vacuum, thus allowing ample time to locate the area of interest and generate images.



The pictured peltier cooling stage device is optional



Accelerating Voltage: 15 kV, UVD image
Sample: Asparagus
Temperature: -20 °C, Vacuum: 100 Pa,
Observing after 10 minutes evacuation



Accelerating Voltage: 15 kV, UVD image
Sample: Asparagus
Temperature: Ambient, Vacuum: 100 Pa,
Observing after 3 minutes evacuation